



germanengineering

In Deutschland: Vertrieb ausschließlich über medizinischen Fachhandel

BM 40 Blutdruckmessgerät

- Vollautomatische Blutdruck- und Pulsmessung am Oberarm
- XL-Display: 83 x 54 mm, Zifferngröße SYS/DIA 21 mm
- Durchschnittswert aller gespeicherten Messwerte und Morgen- und Abendblutdruck der letzten 7 Tage
- 2 x 60 Speicherplätze
- Mit WHO*-Indikator (Einstufung der Messergebnisse gemäß den WHO*-Richtlinien)
- Arrhythmie-Erkennung: erkennt und warnt bei eventuellen Herzrhythmusstörungen
- Medizinprodukt
- Für Oberarmumfänge von 22 - 35 cm
- Datum und Uhrzeit / Abschaltautomatik
- Meldung bei Anwendungsfehler
- Automatische Druckvorwahl und Luftablass
- Mit Blutdruckpass
- Batteriewechselanzeige
- Praktische Aufbewahrungstasche
- Inkl. 4 x 1,5V Batterien AA
- Produktmaße (L/B/H): 119 x 109 x 60 mm
- Produktgewicht: ca. 257 g
- HMV-Nr.: 21.28.01.2100
- XL Manschette (30 - 42 cm) als Zubehör erhältlich Art.-Nr.:162.973
- Netzteil separat erhältlich Art.-Nr.: 071.19
- 3 Jahre Garantie

VE: 4 / Transportkarton: 3 x 4
 EAN-Nr.: 4211125 65815 1
 Art.-Nr.: 658.15

*WHO = World Health Organization (Weltgesundheitsorganisation)

BM 40 Blood pressure monitor

- Fully automatic blood pressure and pulse measurement on the upper arm
- XL display: 83 x 54 mm, digit size SYS/DIA 21 mm
- Average value of all stored measurements and morning and evening blood pressure for the last 7 days
- 2 x 60 memory spaces
- With WHO* indicator (classification of values according to WHO* guidelines)
- Arrhythmia detection: warning in case of possible heart rhythm disturbance
- Medical product
- For upper arm circumferences from 22 - 35 cm
- Date and time / automatic switch-off
- Signal in case of application error
- Automatic pressure preselection and deflation
- With blood pressure passport
- Battery changing indicator
- Practical storage pouch
- Incl. 4 x 1.5V AA batteries
- Product measurements (L/W/H): 119 x 109 x 60 mm
- Product weight: approx. 257 g
- XL cuff (30 - 42 cm) available as an accessory Item no.:162.973
- Mains part available separately Item no.: 071.19
- 3 year guarantee

Sales unit: 4 / Shipping carton: 3 x 4
 EAN no.: 4211125 65815 1
 Item no.: 658.15